

# ASSISTIVE TECHNOLOGY PRODUCT CATALOG 2019



## QUHA ZONO - GYROSCOPIC WIRELESS MOUSE



Quha Zono is the original gyroscoptic mouse designed for special needs.

Quha Zono makes it possible to use a computer when you cannot use an ordinary mouse with your hands. It translates natural head movements into mouse pointer movements – the mouse pointer follows your head movements.

It is an extremely light device, and since it is completely wireless there are endless ways to wear it.

WITH QUHA ZONO YOU CAN DO EVERYTHING YOU CAN DO WITH THE MOST ACCURATE DESKTOP MICE.

IS QUHA ZONO FOR ME?

It provides accurate and intuitive computer access with just small head movements!

CAN I USE QUHA ZONO WITH MY DEVICES?

It is compatible with all available computers and almost all tablets and smartphones! \*)

HOW DO I START USING QUHA ZONO?

There is absolutely zero setup time: open the box, plug it in and turn it on. Instant access to your device!

\*) Apple iOS devices (iPhone & iPad) do not support mice as input devices.



//

I have been waiting for this mouse for 10 years.

Paula, Zono user from Finland

## HEADBAND

Headband is an easy way to wear Quha Zono. It is light and fits comfortably for hours. Headband is adjustable to make it suitable for most users.

Headband is an ideal alternative when using Zono while resting the head on a pillow or a headrest.

Quha Zono is easily attached directly to the Headband, making it very light to wear. The Quha Zono mouse can be charged while attached to the Headband.



## EYEWEAR KIT

Eyewear kit is a spectacle type way to wear Quha Zono. It includes a light and firm eyeglass frame and the Light Eyeglass Clip.



## UNIVERSAL WEAR KIT

Universal Wear Kit enables Quha Zono usage by hand or by foot. It includes a Velcro band and a soft holder for Quha Zono.



## HEAD MOUNT KIT

Head Mount Kit is a smart way to wear Quha Zono. It fits comfortably on most users' heads and can be worn for hours. Putting on and removing Quha Zono is easy, comfortable and fast!

Head Mount Kit includes Light eyeglass clip and Neckband.



## LIGHT EYEGLASS CLIP

Light Eyeglass clip is an ultra light clip to attach Quha Zono safely to your own eyeglasses.



## BASEBALL CAP

Quha Baseball cap includes a classic cap and a white soft holder for Quha Zono.



## QUHA NEMO - SELF OPERATED HEADBAND WITH CHARGING DOCK



NEW!

Quha Nemo is a modified version of the Headband. It includes a charging dock with a possibility to pick it up independently without any assistance.

Nemo is designed so that the user can pick up the Headband with their head and start using it immediately. There is no need for assistance, and after using it, the user can return it to the charging dock where it will automatically start charging.

There are two variants of the Nemo, one with a Sento puff switch and one without. Quha Nemo comes with a Zono mouse, Nemo Headband with or without Sento and the charging dock and charger.

### INDEPENDENT WAY TO WEAR AND USE THE ZONO MOUSE.



#### WHO IS IT FOR?

Nemo is for those who wish to use the Zono mouse independently without any assistance.

#### HOW DO I INSTALL THE NEMO?

Installation is very easy. All you need is a mount with a camera thread. Attach the Nemo to it, plug in the charger and you're good to go.


#### WHAT DOES IT ENABLE?

You can take full control of your computer usage. No assistance required, easy charging and the ability to make mouse clicks. \*)

\*) There are two variants of Nemo. One with Sento puff switch and one without.

Quha Nemo makes it possible to use Quha Zono independently.

## QUHA SENTO - CONTACTLESS PUFF SWITCH



Quha Sento is the world's first contactless puff switch. It is used to perform a short mouse click or dragging together with the Quha Zono gyrosopic mouse.

Mouse clicks are made easily and hygienically with a small puff to the microphone type of mouthpiece.

WITH QUHA SENTO YOU CAN PERFORM  
MOUSE CLICKS EASILY, HYGIENICALLY AND  
WITHOUT HANDS.

---

## QUHA SENTO MP - MULTIPURPOSE CONTACTLESS PUFF SWITCH

Quha Sento MP is a generic multipurpose contactless puff switch. Quha Sento MP enables users to control any assistive technology system in the world. For example environmental control, eye tracking, nurse calls and other systems can easily be controlled by Quha Sento MP. It replaces any ability switch with the same standard connector.

Quha Sento MP includes a flexible sensor arm, a camera thread, a buzzer and configurable sensitivity.



## QUHA PUFO+ - PUFF SWITCH FOR QUHA ZONO

Quha Pufu+ is a wireless puff switch for Quha Zono gyroscopic mouse, making it even more versatile as a headmouse.

Quha Pufu+ puff switch is quickly connected to Quha Zono — no additional power or configurations are needed. The mouthpiece is easily adjustable to each user and replaceable for hygiene.

A little puff to Pufu's mouthpiece will generate a left mouse click and two short puffs a double-click. The dragging can be easily done with a long puff.

Quha Pufu+ comes with an adjustable Headband or Head Mount Kit.



---

## QUHA KAJO BUTTON - ABILITY SWITCH



Quha Kajo Button is a sensitive and durable Assistive Technology switch designed for everyday use. It can be used in any situation that requires switch access.

Quha Kajo Button has been designed to endure hard use and no compromises have been made with the quality of components. It is mountable using the thread inserts in the bottom.

Quha Kajo Button switches have distinctive auditory and tactile feedback. They are available in two different sizes - 25 mm (1 inch) and 50 mm (2 inches).

---

## QUHA KAJO SOFT BUTTON - PILLOW SWITCH

Quha Kajo Soft Button is a soft switch with a large soft surface. Because of the soft foam layer under the cover the Soft Button can be safely pressed with head movement.

The safety pin in the bottom of the Soft Button ensures that the button stays where it is supposed to stay during the night. It can be attached to a pillow or other fabrics. The safety pin is detachable.

The cover of the Soft Button is dirt resistant and skin friendly. The cover can be washed (up to 60°C/140°F) and replaced. Each Soft Button includes one cover.





## QUHA VIDO - EXTERNAL VIRTUAL KEYBOARD AND MOUSE

Quha Vido is an external virtual keyboard and mouse. Quha Vido allows you to use your computer more efficiently, since it frees up your screen space as there is no need to display a virtual keyboard and mouse on your computer screen. Using Quha Vido makes writing emails, browsing the internet and working with your computer faster and more efficient.

The easy and fast setup lets you decide where you want to place the Vido display. You can use Quha Vido together with any wireless computer mouse provided by other manufacturers, although using it in conjunction with the Quha Zono wireless computer mouse gives you the best user experience. With the help of its adjustable stand, you can place the Vido display exactly how you want it and adjust the angle to your liking.

### EFFORTLESS DISPLAY FOR UNLIMITED COMPUTER ACCESS

- Plug-and-play.
- Automatic one-click installation for virtual keyboard and mouse.
- Place the screen where it best suits you.

### MAKES IT POSSIBLE TO USE COMPUTER EVEN FASTER

- Make room for important things on screen by storing keyboard on an external screen.
- External USB monitor for virtual keyboard and mouse.
- Powered from a USB cable.



//

I love the Quha Vido. I can use my computer much faster now

Vido user from Finland



Designed in Finland by Quha

[www.quha.com](http://www.quha.com)

Copyright 2018 Quha oy.

Quha, the Q-mark, Quha Zono, Quha Sento, Quha Pufu+ and Quha Vido are trademarks or registered trademarks of Quha oy.